
CREATE

BALANCE BODY BIO

MULTIFUNCTION DIGITAL SCALE
BÁSCULA DIGITAL MULTIFUNCIÓN

USER MANUAL

[illegible]

[illegible]

INDEX

ENGLISH

Security instructions	6
Buttons	7
Operation	7
Standard references	8
Body fat percentage standard	8
BMI index standard	8
Body hydration percentage standard	8
Body Muscle Percentage And Bone Mass Percentage Standard	8
Basal metabolic rate standard	9

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ENGLISH

Thank you for choosing our scale. Before using the appliance, and to ensure the best use, carefully read these instructions.

The safety precautions enclosed herein reduce the risk of death, injury and electrical shock when correctly adhered to. Keep the manual in a safe place for future reference, along with the completed warranty card, purchase receipt and package. If applicable, pass these instructions on to the next owner of the appliance. Always follow basic safety precautions and accident-prevention measures when using an electrical appliance. We assume no liability for customer failing to comply with these requirements.

SECURITY INSTRUCTIONS

Use the device only as described in the instructions. The company will not be held responsible for any damages which may derive from an improper use.

- Please take off shoes and socks and keep your bare feet in touch with the electrodes when taking body fat %, water %, muscle %, bone substance %, BMI and Calorie measurement.
- Please check the batteries if the scale malfunctions. Put new batteries if necessary.
- Please use soft tissue with alcohol or glass cleaner to clean the surface if it is dirty. No soap or other chemicals. Keep it away from water, heat and extreme coldness.
- The scale is a high precision measuring device. Never jump or stomp on the scale or disassemble it and please handle it carefully to avoid breaking when moving it.
- The scale is for family use only and not suitable for professional use. Body fat %, water %, muscle %, bone substance %, BMI and Calories measured are for reference only. You should consult a doctor when you undertake any diet or exercise program.

Caution: Slippery when wet ! The glass platform may be pretty slippery when wet, thus please make sure both the glass and your feet are dry before use. Never step on the scale platform with wet feet.

- Persons with pacemakers should not use this product.
- Persons with metal implants will be unable to receive accurate readings.
- To ensure accuracy, we suggest you take measurement at the same time every morning.
- Please exchange the current batteries for new ones when “Lo” appears on the LCD screen, indicating the current batteries are out of power.
- “Err” showed on LCD indicates overload.
- The LCD will display “Err1” if a body fat percentage is less than 5% or over 50%.

- This device can be used by children 8 years of age and older, and adults, with reduced physical, sensory or mental capabilities or lack of experience and knowledge, if they have been given supervision or instruction regarding the safe use of the device and understand the hazards involved.

BUTTONS

- SET key: Power on, menu select and confirm
- UNIT key: Unit change (kg /lb/ st:lb)
- Up key: To move up the menu
Down key: To move down the menu

OPERATION

- Batteries are included in the box, insert them into the battery compartment on the back of the scale.
- Please use the scale on a hard, flat and clean surface. Weight units can be adjusted and chosen from Kg, Lb, and ST: LB from the button on the back of the scale.

For weight only

1. Press the base and it will get ready for use when the screen shows “0.0kg” after “88.8” and you can press Up and Down buttons to switch within different units (kg/lb/st:lb).
2. Step on the scale, position your feet evenly near the center and keep your body still until the scale stabilizes your weight. The scale will turn off automatically if there is no further use.

For weight, body fat %, water %, muscle %, bone substance %, BMI and Calories.

1. Long press SET key to set your personal data when the screen shows 0.0 kg. Press Up or Down keys to select a number from 0 to 9 as your personal code when the code number flashes on the screen and once you reach your target number, please press SET again to confirm.
 2. Press Up or Down buttons to choose male or female profile when the gender symbol flashes on the screen and then press SET to confirm.
 3. Press Up or Down buttons to reach your target height (100-240 cm) when the number and cm icon flash on the screen and then press SET to confirm.
 4. Press Up or Down buttons to reach your target age (10-80) when the word age and the number flash on the screen and then press SET to confirm.
 5. Step on the scale gently and keep your bare feet in touch with the electrodes. Stand still and maintain an upright posture to get your weight, body fat percentage, water percentage, bone substance percentage, BMI and total required calories intake readings. All values will be displayed for a few seconds, repeated twice, and the scale will be turned off automatically following no further usage.
- The personal data will be stored once it is set. Confirm your user number everytime you use the scale, your information will be shown automatically, there is no need to set it up again.

STANDARD REFERENCES

- These charts should be used **ONLY FOR REFERENCE**.
- Consult your doctor before starting any diet or exercise program.

BODY FAT PERCENTAGE STANDARD

Rating	Standard for men					Standard for women				
Age	20-29	30-39	40-49	50-59	60+	20-29	30-39	40-49	50-59	60+
Low	<13	<14	<16	<17	<18	<19	<20	<22	<23	<24
Normal	14-20	15-21	17-23	18-24	19-25	20-26	21-27	23-29	24-30	25-31
Moderately high	21-23	22-24	24-26	25-27	26-28	27-29	26-28	30-32	31-33	32-34
High	>23	>24	>26	>27	>28	>29	>30	>32	>33	>34

BMI INDEX STANDARD

Underweight	Normal weight	Slightly overweight	Overweight	Severe overweight	Extremely overweight
<18.5	18.5-24.9	25.0-29.9	30.0-34.9	35.0-39.9	>40

BODY HYDRATION PERCENTAGE STANDARD

Gender	Age	Slight	Low	Standard	High
Male	10-20	<53	53-57	57-70	>70
	20-40	<50	50-54	54-62	>62
	40-60	<48	48-53	53-59	>59
	60-80	<44	44-49	49-57.5	>57.5
Female	10-20	<52	52-55	55-63	>63
	20-40	<49	49-52	52-60	>60
	40-60	<45	45-49	49-56	>56
	60-80	<39	39-43	43-52.5	>52.5

BODY MUSCLE PERCENTAGE AND BONE MASS PERCENTAGE STANDARD

Muscle percentage				
Gender	Slight	Low	Standard	High
Male	<35	35-39	39-49	>49
Female	<30	30-35	35-44	>44
Bone mass percentage				
Gender	Slight	Low	Standard	High
Male	<4.2	4.2-4.5	4.5-5.1	>5.1
Female	<3.6	3.6-3.9	3.9-4.3	>4.3

BASAL METABOLIC RATE STANDARD

Age	Female	Male
1-2	700 (kcal/day)	700 (kcal/day)
3-5	860 (kcal/day)	900 (kcal/day)
6-8	1000 (kcal/day)	1090 (kcal/day)
9-11	1180 (kcal/day)	1290 (kcal/day)
12-14	1340 (kcal/day)	1480 (kcal/day)
15-17	1300 (kcal/day)	1610 (kcal/day)
18-29	1210 (kcal/day)	1550 (kcal/day)
30-49	1170 (kcal/day)	1500 (kcal/day)
50-69	1110 (kcal/day)	1350 (kcal/day)
70-	1010 (kcal/day)	1220 (kcal/day)

In compliance with Directives: 2012/19/EU and 2015/863/EU on the restriction of the use of dangerous substances in electric and electronic equipment as well as their waste disposal. The symbol with the crossed dustbin shown on the package indicates that the product at the end of its service life shall be collected as separate waste. Therefore, any products that have reached the end of their useful life must be given to waste disposal centres specialising in separate collection of waste electrical and electronic equipment, or given back to the retailer at the time of purchasing new similar equipment, on a one for one basis. The adequate separate collection for the subsequent start-up of the equipment sent to be recycled, treated and disposed of in an environmentally compatible way contributes to preventing possible negative effects on the environment and health and optimises the recycling and reuse of components making up the apparatus. Abusive disposal of the product by the user involves application of the administrative sanctions according to the laws.

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