
# CREATE

## Balance Body Bamboo

---

User manual | Manual de instrucciones

# CREATE

## INDEX

---

### ENGLISH

Safety warnings	6
Operating instructions	7
Automatic switch off	7
Overload indicator	7
Battery change indicator	7
Cleaning and maintenance	7

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Thank you for choosing our scale. Before using the appliance, and to guarantee the best use, read these instructions carefully.

The safety precautions listed here reduce the risk of fire, electric shock, and injury when followed correctly. Please keep the manual in a safe place for future reference, as well as the warranty, sales receipt and box. If applicable, give these instructions to the future owner of the appliance. Always follow basic safety instructions and risk prevention measures when using an electrical appliance. The manufacturer will not be held responsible for any damage resulting from the failure of the user to follow these instructions.

## SAFETY WARNINGS

---

Use the device only as described in the instructions. The company will not be held responsible for any damages which may derive from an improper use.

- Please check the batteries if the scale malfunctions. Replace them if necessary.
- Please use soft tissue with alcohol or glass cleaner to clean the surface if it is dirty. No soap or other chemicals. Keep it away from water, heat and extreme coldness.
- The scale is a high precision measuring device. Never jump or stomp on the scale or disassemble it and please handle it carefully to avoid breaking when moving it.
- The scale is for family use only and not suitable for professional use. You should consult a doctor when you undertake any diet or exercise program.

**Caution:** Slippery when wet! The scale platform may be pretty slippery when wet, thus please make sure both the scale platform and your feet are dry before use. Never step on the scale platform with wet feet.

- Please use the scale on a hard and flat surface. Do not use on a carpet or soft surface.
- Carefully step on the platform of the scale. Stand as stationary as possible until the weight reading is shown and locked on the display.
- The scale will be automatically off if there is no operation for a while.
- This device can be used by children 8 years of age and older, and adults, with reduced physical, sensory or mental capabilities or lack of experience and knowledge, if they have been given supervision or instruction regarding the safe use of the device and understand the hazards involved.

## OPERATING INSTRUCTIONS

---

1. Press the platform centre and remove your foot.
2. "0.0" will be displayed.
3. The scale will switch off and now it is ready for use.
  - This initialisation process must be repeated if the scale is moved. At all other times step straight on the scale.
4. Gently step on the scale. Put your feet evenly on the scale without moving and wait whilst the scale computes your weight. Your correct weight will then be displayed in the following way: 80.3 kg, for example.

## AUTOMATIC SWITCH OFF

---

- When you step off the scale, it will automatically switch off. (If no weighting is carried out while the zero is being displayed, the scale will automatically switch off after approximately 6 seconds).

## OVERLOAD INDICATOR

---

- When the scale is over the maximum capacity, the following will appear on the screen: ERR

## BATTERY CHANGE INDICATOR

---

- When the batteries need replacing, the following will appear on the screen: LO
- Please open the back board of the scale, then take the batteries carefully out of the battery compartment. Please push the metal sheet outwards by your fingernails or a sharp instrument.
- The battery will pop up. Do not bend the contacts back while you do this, or they may break. Dispose of the old batteries in a safe place.

## CLEANING AND MAINTENANCE

---

- Placing your scale on a hard flat surface will ensure the greatest accuracy and repeatability.
- Always weigh yourself on the same scale placed on the same floor surface. Uneven floors may affect the readings.
- Always weigh yourself without clothing or footwear, and before meals. First thing in the morning is the best time.
- Clean your scale with a damp cloth but do not allow water to get inside the scale. Do not use any chemical cleaning agents.
- Do not allow your scale to become saturated with water as this can damage the electronics components.
- Do not store the scale in an upright position when not in use as this will drain the batteries.
- Remove the batteries if the scales are not going to be used for a long period of time.
- Treat your scale with care as it is a precision instrument, do not drop it or jump up or down on it.

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In compliance with Directives: 2012/19/EU and 2015/863/EU on the restriction of the use of dangerous substances in electric and electronic equipment as well as their waste disposal. The symbol with the crossed dustbin shown on the package indicates that the product at the end of its service life shall be collected as separate waste. Therefore, any products that have reached the end of their useful life must be given to waste disposal centres specialising in separate collection of waste electrical and electronic equipment, or given back to the retailer at the time of purchasing new similar equipment, on a one for one basis. The adequate separate collection for the subsequent start-up of the equipment sent to be recycled, treated and disposed of in an environmentally compatible way contributes to preventing possible negative effects on the environment and health and optimises the recycling and reuse of components making up the apparatus. Abusive disposal of the product by the user involves application of the administrative sanctions according to the laws.

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